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RUS

Dear customer!

We congratulate you on purchasing home appliance of „MAGNIT“ trademark. We are sure that this unit will correspond to the highest requests for the quality and will become the real reliable support in your home cares. Before starting using this appliance, please carefully read the instruction manual and follow the rules of using. In this case the trademark „MAGNIT“ will assure the long life of the equipment you have bought.

Thank you for choosing us!

## Description

Note: before starting using this appliance please carefully read the instruction manual.

### Safety precautions

1. Congratulation for having procured this body fat scale as monitoring device of your body fat %, as well as body hydration %. The principle of this device is based on the measurement of body's electronic impedance known as 'bio-electronic Impedance Analysis' (BIA), a factor related to the body fat & hydration % and in concurrence with other physical data(Age, Gender, Body Height).
2. The evaluated readings of body fat & hydration % are only reference value. It is not recommended for use of any medical purpose. Please always consult your doctor for advice of medication or diet in order to achieve ideal weight level.
3. Pacemaker users/wearers, or other medical implants with electronic components, are not recommended to use this scale. Implant performance may be jeopardized and impaired.
4. To prevent interference, do not use the device in environment near equipment with strong electromagnetic field. Always keep a good distance from such equipments when operating the scale, or simply turn off those equipments radiating strong electro-magnetic fields before operating.


## FEATURES / SPECIFICATIONS:

### 1. Memory, Gender, Height, Age setti.

PARAMETER	DEFAULT SETTING	RANGE
Memory	1	1 to 8
Gender	Male	Male/Female
Height	165 cm	100 cm - 250 cm
Age	25	10-99

2. Max. Weighing Capacity: 150 kg/330 lb/24 stone
3. Min. Weighing at: 2,0 kg/4,4 lb
4. Weight Measuring Division: 0,1 kg/0,2 lb
5. Kg/Lb/Stone convertible
6. Body Fat measure range: 4,0-45,0%
7. Body Fat differential resolution in %: 0,1%
8. Hydration measure range: 37,8-66,0%
9. Hydration differential resolution in %: 0,1%
10. Main LCD digit size: 1.0"(23mm)
11. Auto-zero, auto-off
12. Low battery power indication ('LO')/Overload indication ('Err')
13. Power: 1x9 V alkaline battery.

### LCD SEGMENTUM / ICONS ILLUSTRATIONS:

### INITIAL PREPARATION:

1. Install 1 x 9 V battery into the compartment underneath the main scale unit.
2. Place the unit on a dry, flat and solid ground surface, ensure the unit is standing stable. Avoid carpet, or wavy ground surface.

ON/SET = : On/Set switch

START = : Start switch


: UP arrow press switch

: Down arrow press switch

## FAT MONITOR SCALE OPERATION:

### Initial Set Up:

1. Press button 'on/set' Key to wake up the unit. All LCD segmentum light up briefly and prompt to setting mode starting from Memory selection, By pressing button 'on/set' key again will take you through cycle of setting all parameters. See following illustration:

- a. Press Up / Down Arrow (▲▼) button to advance reverse to the memory station which you would like to choose or has been designated to you. Press (on/set) button to confirm setting, LCD will prompt to next parameter setting (Gender) at the same time.
  - b. Press Up / Down Arrow (▲▼) button to advance / reverse to choose the gender of your type. Press (on/set) button to confirm setting. LCD will prompt to next parameter setting (Height) at the same time.
  - c. Press Up / Down Arrow (▲▼) button to advance / reverse to the adjust to your height level. Press (on/set) button to confirm setting. LCD will prompt to next parameter (Age) at the same time.
  - d. Press Up / Down Arrow (▲▼) button to advance/reverse to adjust to your age. Press (on/set) button to confirm setting. LCD will prompt to Memory setting. All setting complete!!!  
At this point, you can continue to press (on/set) button a few consecutive times to re-cycle through the setting mode to review and assure all the parameters are correctly set.
  - e. Press (start) button to exit setting mode and enter normal Fat monitor operating mode. LCD will prompt to '0.0 kg' (or '0.0' lb). You are ready to weigh yourself, and let the scale compute your Fat %, Hydration %.
2. Step on the scale gently (ensure to remove your socks & shoes, wipe clean your feet & the weighing platform to get the best conductive effect). While stepping on, align your feet to the metal conductive pads, two for each foot. Ensure your feet are standing apart not touching each other (to avoid shortcircuit).  
If you step on the scale too late after the LCD auto-off, simply press the (on/set) button again, wait till LCD prompt to show your current settings. Then press (start) button again to reenter normal fat monitor scale operation mode.

3. As soon as you step on the scale, the upper side area of the LCD will immediately display incrementing digits till it stabilises, locked & flashes, indicating your current weight.
4. Within a few seconds, the scale will automatically compute your Fat % and Hydration %. In the lower left corner, LCD will prompt to show the measuring process and result of your Fat % (icon "BF") and hydration% (icon "BW"). Subsequently, LCD will prompt back and forth alternately between value of Fat % and hydration % for 5 times, and then automatically shut down.

#### **NOTE:**

Whenever pressing **(on/set)** button, turning on the unit for fresh operation, you can always go through the same steps from **Illustration (a)** to **(e)** to verify or adjust all personal data and previous measured records, before measuring for new records.

#### **NORMAL WEIGHING SCALE OPERATION:**

**To use as General Weighing Scale (Disregard of any parameter setting for Fat / hydration%).**

1. Turn on the unit by slight tap on the lower right corner of the platform of the scale. LCD will wake up and prompt to General Weighing Mode, showing "0.0 kg" (or "0.0 lb")
2. Step on the scale gently (Do not jump on the scale). As soon as you step on the scale, LCD will immediately display incrementing digits till it stabilises followed by flashes and then locked digits, showing your weight value. After a few seconds, LCD will automatically shut down.

#### **ERROR & OTHER LCD MESSAGE:**

LCD may show occasionally error message 'Err'. Simply turn off the unit for a brief moment, or remove and re-install battery as a total re-set. Then turn on again to resume operation.

1. All segmentum lit up unanimously:  
Indicating all segmentum performing correctly.
2. "-——": Strain gauge sensor zero point setting too low below lower range. Can not operate properly.
3. "Err ": Overload, or strain gauge sensor zero point too high beyond upper range.
4. "LO": Main unit battery power running low.
5. "Err-h": The computed **Fat** % rate too high above range of upper limit.
6. "Err-L": The computed **Fat** % rate too low below range of lower limit.

## MAINTAINING & CLEANING PROPER USAGE:

1. Use a dry, soft fabric to clean the device. Or take a lightly damp soft fabric with mild detergent to wipe clean the device. Do not use strong chemical solution as alcohol, benzene for cleaning.
2. It could be slippery when the unit platform surface is wet, or semi-damp. Ensure the weighing platform surface is dry before using for safety reason.
3. Always weigh/use the same scale on a hard flat surface. Do not use on carpets, mats. Remove the battery if the scale is not going to be used for a long period of time.
4. Do not drop the device, or jump on the device when using. Avoid strong impact and sudden shake as it is a precision measuring device. The construction is delicate and sensitive to all unwanted external impact.
5. Do not insert any exterior objects into any openings on the unit.
6. Do not attempt to open the unit to fix & repair. This will invalidate the guarantee.
7. If fail to turn on the unit, check whether battery had been installed, or battery power is too low. Insert or replace new battery. Or contact your dealer for consultation.
8. If there are some error on the display or can not turn off for a long time, please take off the battery for about 3 seconds. Then encase it again to remove the soft-fault.
9. Improper use of the scale unit for purpose other than normal weighing, or exposing the unit to electro-magnetic wave which leads to loss of data will not be covered by our warranty/guarantee, nor our company be liable for any of such loss or mishap.
10. Our company will not be responsible for mis-use of the scale unit for other purpose that incurs physical damages to the user or any other third party user.

## REFERENCE INFORMATION:

The evaluated body fat result is based on a standardised measuring method. It is calculated by a correlated formula developed on basis of extensive scientific studies. The hydration contact in the lean muscle affects measuring result. Individual habits and lifestyle will affect balance of hydration level. Other factors, like following may affect hydration balance:

- a. After a bath, readings can be lower.
- b. After a meal or drinking readings can be higher.
- c. During ladies menstruation period can be fluctuated/inconsistent.
- d. Loss of hydration by illness, strenuous physical activities that cause heavy sweats and dehydration.
- e. Table below gives estimated physical status relating to the indicated **Fat % Hydration%** from the scale:

Age	Fat% Hydration%		Fat% Hydration%		Status	
	Female	Female	Male	Male		
below 30	4.0-16.0%	66.0-57.8%	4.0-11.0%	66.0-61.2%	too lean	
	16.1-20.5%	57.7-54.7%	11.1-15.5%	61.1-58.1%	lean	
	20.6-25.0%	54.6-51.6%	15.6-20.0%	58.0-55.0%	normal	
	25.1-30.5%	51.5-47.8%	20.1-24.5%	54.9-51.9%	high very	
	30.6-45.0%	47.7-37.8%	24.6-45.0%	51.8-37.8%	high	
above 30	4.0-20.0%	66.0-55.0%	4.0-15.0%	66.0-58.4%	too lean	
	20.1-25.0%	54.9-51.6%	15.1-19.5%	58.3-55.3%	lean	
	25.1-30.0%	51.5-48.1%	19.6-24.0%	55.2-52.3%	normal	
	30.1-35.0%	48.0-44.7%	24.1-28.5%	52.2-49.2%	high very	
	35.1-45.0%	44.6-37.8%	28.6-45.0%	49.1-37.8%	high	

**PLEASE READ FOLLOWING NOTES CAREFULLY FOR PROPER UNDERSTANDING OF THE APPLICATION/SIGNIFICANCE OF THE DEVICE:**

1. The **Fat %/Hydration%** content evaluated by the Fat Monitor Scale is only for reference. The physical outlook of body does not necessarily gives congruent **Fat %/Hydration%** readings, nor the same physically built persons of same age and sex should give the same readings, Please consult one's doctor to determine the best condition of body **Fat/Hydration** percentage that one should achieve.
2. Body Fat Scale is not suitable for people of following conditions to use :
  - \* Children under 10, and adults over 99
  - \* People with symptoms of fever, edema, osteoporosis.
  - \* People on dialysis treatment.
  - \* People taking cardiovascular medicine.
  - \* Being Pregnant
  - \* Athletes who undertake more than 10 hours of intensive sport activity.
  - \* Those with heart pulse rate under 60 beats/minute.

**Technical characteristics**

Weighing diapason: 0 - 150 kg  
 Accuracy: 0.1 kg  
 Power supply : one 9-voltaic battery power (does not enter into a delivery package), without a power supply network.

The life-time of product is 2 years.  
 You can find full information about the appliances of "Magnit" trademark on [www.magnit.net](http://www.magnit.net).

The manufacturer reserves the right to bring modifications to the design, construction and completeness of the product without prior notice.

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